

swash

another Arion brand

Swash® SPECIALS hair-care

'Washing without water', the innovative and user-friendly washing concept



Problem

Maintaining levels of quality and consistency of care given by carers are critical factors in the washing of patients' hair. In fact, this is impossible when using the traditional washing method. The washing process can be adversely affected by a number of different factors, such as the temperature of the water, the amount of water, the quality of the shampoo, the quantity of shampoo, etc. Furthermore, there is static load, a major source of physical strain when caregivers have to work for long periods in difficult postures. Although static load is also common outside basic care situations, it is particularly prevalent when carers have to wash, shower and bandage patients.

This problem is exacerbated by the time pressure caregivers have to work under. This means there is less time for other important indicative and supervisory tasks associated with the washing of patients' hair.

Solution

Swash SPECIALS hair-care guarantees high levels of quality and consistency in care by perfectly combining its various components: the ingredients of the lotion, the amount of lotion and the conditioner. The lotion has been proven to have a positive effect on skin irritations. The risk of cross-contamination, wastage and spills is also reduced on account of the fact that Swash SPECIALS hair-care is used on a patient-by-patient basis.

The hair-care lessens the physical strain on the caregiver by reducing the number of washing actions to 1. A traditional hair wash, including preparation, can take anything up to 25 minutes, whereas Swash SPECIALS hair care can cut this down to just 8 minutes.

From an ergonomic and cost-saving point of view – in addition to the high levels of consistency and quality achieved – the 'washing without water' concept offered by the Swash SPECIALS hair care is the best there is, certainly for bed-ridden patients, terminal patients or patients who are treated with cytostatics.

Target group

Any carer who has to wash the hair of a patient.

Swash® SPECIALS hair-care

Top quality in care

Product benefits

- Consistency and high levels of quality in care
- Effortless but thorough washing of bed-ridden patients' hair
- Reduces physical strain on carers and patients alike
- Is less time-consuming than the traditional washing method (approx. 75% reduction)
- Reduces the risks of cross-contamination due to less contact
- Skin and eye irritation test available
- Can have a positive effect on absenteeism
- PH-neutral, alcohol and latex free

