

swash

another Arion brand

Swash® PLATINUM wipes

Top quality in care / 'Washing without water', the innovative and user-friendly washing concept



Problem

Maintaining levels of quality and consistency of care given by carers are critical factors in the washing of patients. In fact this is impossible when using the traditional washing method. The washing process can be adversely affected by a number of different factors, such as the temperature of the water, the amount of water, the quality of the soap, the quantity of soap, etc.

Furthermore, there is static load, a major source of physical strain when carers have to work for long periods in difficult postures. Although static load is also common outside basic care situations, it is particularly prevalent when caregivers have to wash, shower and bandage patients.

Other serious problems faced by caregivers when washing patients are skin irritations and infections which are only exacerbated when the work is carried out under the pressure of time. This means there is less time for other important indicative and supervisory tasks associated with the washing of patients.

Solution

Swash PLATINUM wipes guarantee high levels of quality and consistency perfectly combining its various components: the composition of the lotion, the amount of lotion and the packaging. The lotion has a positive effect on skin irritations and infections. Furthermore, there is a reduced risk of cross-contamination due to the fact that Swash PLATINUM wipes are used on patients individually.

The wipes lessen the physical strain on the carer by reducing the number actions by 4 to 1. The duration of the wash is also curtailed.

A traditional bed bath, can take anything up to 25 minutes, whereas Swash PLATINUM this down to just 12 minutes.

From an ergonomic and cost-saving point of view – in addition to the consistency and quality achieved – the 'washing without water' concept Swash PLATINUM wipes is the best there is, certainly for bed-ridden patients or patients who are treated with cytostatics.

Target group

Any carer who washes patients

Swash® PLATINUM wipes

Top quality in care

Product benefits

- Consistency and high levels of quality in care
- Effortless but thorough washing of bed-ridden patients
- Reduces physical strain on carers and patients alike
- Is less time-consuming than the traditional washing method (approx. 50% reduction)
- Reduces the risks of cross-contamination due to less contact
- Skin and eye irritation test available
- Can have a positive effect on absenteeism
- PH-neutral, alcohol and latex free
- Aid in ADL training of patients
- Re-sealable packaging for optimum use
- Stays warm for longer than water in traditional washing



AEB10-009 (PB Swash PLATINUM wipes Export eng)